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Dátum poslednej revízie: 2. 1. 2024

SUMMARY AND EXPLANATION

PYR reagent is used for generation of colour in test PYR in diagnostic kits. Detailed description of work with reagent PYR and interpretation of result of test is part of particular kit manual.

PRINCIPLE

Substrate L – pyrrolidonyl – b – naftylamid is transformed by acting of L – pyrrolidonylmarylamidase to colourless L – pyrrolidon and b - naftylamin. Free b – naftylamin is visualized by PYR reagent by creation of red coloured Shiff base.

DIAGNOSTIC KIT CONTAINS

- PYR reagent min. 10 ml (Ref. 3003)
- Dropper, cap
- information leaflet

WARNINGS AND SAFETY PRECAUTIONS

- **For in vitro diagnostics use.**
- **For professional use only**
- Follow the instructions exactly!
- All samples and inoculated products have to be considered as potentially infectious and obvious safety regulations of your country must be respected while handling.
- Do not use expired product.
- Before use, check if the packaging is intact. Do not use damaged kit.
- Leftovers of expired or devaluated reagent liquidate according to regulations of your country. Packaging materials are recyclable.

STORAGE

PYR reagent can be stored in dark in temperatures between +2 and +8°C. Expiration date is showed on label of each package.

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