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Dátum poslednej revízie: 20. 5. 2022

SUMMARY AND EXPLANATION

HIP reagent is used for generation of colour in test HIP in diagnostic kit. Detailed description of work with reagent HIP and interpretation of result of test is part of particular kit manual.

PRINCIPLE

Substrate sodium hippurate is hydrolysed to end products, glycine and benzoic acid, which is detected by HIP reagent.

DIAGNOSTIC KIT CONTAINS

- HIP a - solution min. 18 ml in glass bottle (Ref. 3006)
- HIP b reagent - solid for reagent in plastic bottle (Ref. 3006)
- dropper
- information leaflet

PREPARING OF REAGENT TO USE

Pour solution "a" to solid "b", then put on dropper. After dissolving of solid, reagent is prepared for use.

WARNINGS AND SECURITY PRECAUTIONS

- **For in vitro diagnostics use.**
- **For professional use only**
- Follow the instructions exactly!
- All samples and inoculated products have to be considered as potentially infectious and obvious safety valid regulations must be respected while handling.
- Before use, check if the packaging is intact. Do not use damaged kit.
- Leftovers of expired or devaluated reagent liquidate according to regulations of your country. Packaging materials are recyclable.

Kit contains mixture of acetone/butanol (50 % / 50 %) and ninhydrin (100 % of solid).

H225 Highly flammable liquid and vapour

H302 Harmful if swallowed

H315 Causes skin irritation

H318 Causes serious eye damage

H335 May cause respiratory irritation

P210 Keep away from heat/sparks/open flames/hot surfaces - No smoking

P242 Use only non-sparking tools.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P270 Do not eat, drink or smoke when using this product.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301+312 IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell

P302+P352 IF ON SKIN: Wash with soap and water

P304+P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

P305+P351+P338 IF IN EYES: Rinse continuously with water for several minutes. Remove contact lenses if present and easy to do continue rinsing

P332+P313 IF skin irritation occurs: Get medical advice/attention.

P370+P378 In case of fire: Use water or sodium bicarbonate for extinction.

P403+P233+P235 Store in a well-ventilated place. Keep container tightly closed. Keep cool

P405 Store locked up.

STORAGE

HIP reagent can be stored in dark at temperatures between +2 and +8°C. Expiration date is showed on label of each package.

Reagent prepared by mixing of components can be stored in dark in temperatures between + 2°C and + 8°C for 5 months or till date of expiration – if expiration date finishes earlier.

Warning: When preparing solution of ninhydrin avoid contact with skin. Contact leads to formation of spots that are difficult to remove.

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H225 Highly flammable liquid and vapour
H302 Harmful if swallowed
H315 Causes skin irritation
H318 Causes serious eye damage
H335 May cause respiratory irritation

P210 Keep away from heat/sparks/open flames/hot surfaces - No smoking

P242 Use only non-sparking tools.

P261 Avoid breathing dust/fume/gas/mist/vapours/spray.

P270 Do not eat, drink or smoke when using this product.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301+312 IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell

P302+P352 IF ON SKIN: Wash with soap and water

P304+P340 IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

P305+P351+P338 IF IN EYES: Rinse continuously with water for several minutes. Remove contact lenses if present and easy to do continue rinsing

P332+P313 IF skin irritation occurs: Get medical advice/attention.

P370+P378 In case of fire: Use water or sodium bicarbonate for extinction.

P403+P233+P235 Store in a well-ventilated place. Keep container tightly closed. Keep cool

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