SUHRN A VYSVETLENIE

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Dátum poslednej revízie: 3. 8. 2017

SUMMARY AND EXPLANATION

VP reagent is auxiliary agent for test VP (Voges – Proskauer test). Reagent shows colour change by creating of acetoín in tests of diagnostic kits and strips produced by DIAGNOSTICS s.r.o.

PRINCIPLE

Glucose is fermented by microorganisms to acetoín. Reagent VP 1 - ethanol solution of α -naftol, deepens final colouring of solution and simultaneously elevates accuracy of reaction. α - naftol catalyses change of acetoín in alkaline environment (VP 2) to diacetyl. KOH reacts with peptone. Condensation of diacetyl with guanid groups of peptone in alkaline environment leads to formation of pink colour, what represents positive result of test.

DIAGNOSTIC KIT CONTAINS

- reagent VP 1a – substrate (Ref. 3004)
- reagent VP 1b – solvent min. 18 ml (Ref. 3004)
- reagent VP 2 min. 18 ml (Ref. 3004)
- dropper
- information leaflet

PREPARING OF REAGENT TO USE

Pour VP 1a to VP 1b. After dissolving of substrate, reagent is ready to use. Mark the label of reagent with the preparation and expiration date. Reagent VP 2 is ready for use.

WARNINGS AND SAFETY PRECAUTIONS

- **For *in vitro* diagnostics use.**
- **For professional use only.**
- Follow the instructions exactly!
- All samples and inoculated products has to be considered as potentially infectious and obvious safety regulations of your country must be respected while handling
- Do not use expired product.
- Before use, check if the packaging is intact. Do not use damaged kit.
- Leftovers of expired or devaluated reagent liquidate according to regulations of your country. Packaging materials are recyclable.

Kit contains ethanol (5 % of volume), potassium hydroxide (5 % of volume), α -naftol (5 % of volume)

H302 Harmful if swallowed
H315 Causes skin irritation
H319 Causes serious eye irritation

P280 Wear protective gloves/protective clothing/eye protection/face protection.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P302+P352 IF ON SKIN: Wash with plenty of soap and water
P332+P313 If skin irritation occurs: Get medical advice/attention.
P337+P313 If eye irritation persists: Get medical advice/attention.

STORAGE

Store VP reagent in dark, in temperatures between +2 and +8°C. Expiration date is showed on label of each package.

Prepared reagent VP1 (VP 1a to VP 1b) store in dark at temperatures between +2 and +8°C.

Use prepared reagent till 5 months from date of preparation, or till date of expiration, in case that expiration date passes earlier.

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