## SK

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## EN

### SUMMARY AND EXPLANATION

TRB is detection strip used for differentiation of bacteria belonging to genus *Moraxella* and *Neisseria*. Detection strip changes colour to yellow on base of metabolism of examined sample.

### PRINCIPLE

Breakdown of substrate by enzymatic hydrolysis is attended by colour change and creation of yellow colouring, this indicates positive result of TRB test.

### DIAGNOSTIC KIT CONTAINS:

- 50 test strips for 50 determinations
- information leaflet

### WARNINGS AND SAFETY PRECAUTIONS

- **For *in vitro* diagnostics use.**
- **For professional use only.**
- Follow the instructions exactly!
- Used strips should be considered as potentially infectious and obvious safety valid regulations must be respected when handling.
- Observe common safety measures according to the regulations of your country.
- Before use, check if the packaging is intact. Do not use damaged kit.

### STORAGE

At dark place at +2 to +8 °C. Expiration date is placed on the label of each package.

### REQUIRED MATERIALS

Tube with 0,5 ml of saline solution

### RECOMMENDED PROCEDURE:

1. Prepare thick bacterial suspension (≥ 2 McF) from fresh 24 hour culture.
2. Put in diagnostic strip and dip zone into suspension.
3. Incubate in thermostat by temperature 35 ± 2°C for 2 - 4 hours.
4. Evaluate result of TRB test.

### INTERPRETATION OF REACTION

**Yellow colouring indicates positive result of TRB test.**

Bacterial strain	TRB
<i>Moraxella catarrhalis</i>	+ (yellow)
<i>Neisseria sp.</i>	- (pink)

For the need of own proof function, use following recommended bacterial strains:

Control strain	CCM / ATCC	TRB
<i>Moraxella catarrhalis</i>	4391 / -	+ (yellow)
clear water	-	- (pink)

CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420549491430, e-mail: ccm@sci.muni.cz

Kit components doesn't contain hazardous substances.

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