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Dátum poslednej revízie: 20. 5. 2022

SUMMARY AND EXPLANATION

Diagnostic discs BAC 10 with content of 10 IU of bacitracin are designed for preliminary identification of haemophilus in primoculture on base of its high resistance to bacitracin with concentration 10 IU, against other microorganisms.

PRINCIPLE

Haemophilus species are able to grow in zone with high concentration of bacitracin (10 IU) by attendant diffusion of grow factors, which are produced by necessarily present staphylococcus.

DIAGNOSTIC KIT CONTAINS

- tube with diagnostic discs - min. 500 pc (Ref. 6023-500), min. 50 pc (Ref. 6023-50)

WARNINGS AND SECURITY PRECAUTIONS

- **For in vitro diagnostics use.**
- **For professional use only**
- Follow the instructions exactly!
- All samples and inoculated products have to be considered as potentially infectious and obvious safety valid regulations must be respected while handling.
- Do not use product out of expiration.
- Before use, check if the packaging is intact. Do not use damaged kit.
- Leftovers of expired or devaluated product liquidate according to regulations of your country. Packaging materials are recyclable.

STORAGE

Keep antibiotic discs in dark at temperatures between +2 and +8 °C

REQUIRED MATERIALS

- Petri dish with blood agar
- tampons, loops, burner, tubes and other basic equipment of microbiological laboratory
- culture of *Staphylococcus aureus*

RECOMMENDED PROCEDURE

Spread tested culture on Petri dish with blood agar. Spread line of culture of *Staphylococcus aureus* to middle and put there sterily disc BAC 10, press it gently, so that whole surface of disc is in contact with medium. Incubate for 18-24 h in CO₂ atmosphere at temperature 37°C.

EVALUATION AND INTERPRETATION

Culture which grows in inhibition zone of disc and near staphylococcus line can be considered as culture of haemophilus.

Functionality of diagnostic discs can be verified by these control strains:

Bacterial strain	CCM / ATCC	Inhibition zone
<i>Haemophilus influenzae</i>	4456 / 49247	grow
<i>Streptococcus pneumoniae</i>	4501 / 49619	> 17 mm

CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420549491430, e-mail: ccm@sci.muni.cz

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