Antibiotický diagnostický disk / Antibiotic diagnostic disc / Antibiotikus diagnosztikai korong

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Dátum poslednej revízie: 20.05.2022

SUMMARY AND EXPLANATION

Diagnostic discs NEO with content of 30 µg neomycin are designed for determination of microorganism to neomycin sensitivity.

PRINCIPLE

Microorganisms are differentiated according to resistance to disc with content of neomycin.

DIAGNOSTIC KIT CONTAINS

- tube with diagnostic discs - min. 500 pcs (Ref. 6040)

WARNINGS AND SECURITY PRECAUTIONS

- **For in vitro diagnostics use.**
- **For professional use only**
- Follow the instructions exactly!
- All samples and inoculated products have to be considered as potentially infectious and obvious safety valid regulations must be respected while handling.
- Do not use product out of expiration.
- Before use, check if the packaging is intact. Do not use damaged kit.
- Leftovers of expired or devaluated product liquidate according to regulations of your country. Packaging materials are recyclable.

STORAGE

Keep antibiotic discs in dark at temperatures between +2 and +8 °C

REQUIRED MATERIALS

- Petri dish with Mueller-Hinton agar
- tampons, loops, burner, tubes and other basic equipment of microbiological laboratory

RECOMMENDED PROCEDURE

Prepare suspension of tested strain in saline (cca 0,5 McF). Spread suspension on whole area of Petri dish with sterile tampon soaked in prepared suspension. Place sterilely disc NEO on inoculated Petri dish, press it gently, so that whole surface of disc is in contact with medium. Incubate for 18-24 h at temperature 37°C.

EVALUATION AND INTERPRETATION

Inhibition zone


≤12 mm	13-16 mm	≥ 17 mm
susceptible	intermediate	resistant

Functionality of diagnostic discs can be verified by these control strains:

Bacterial strain	ATCC / CCM	Inhibition zone
<i>E. coli</i>	25922 / 3954	17-23 mm

CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420549491430, e-mail: ccm@sci.muni.cz

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