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Dátum poslednej revízie: 20. 5. 2022

**SUMMARY AND EXPLANATION**

Diagnostic discs NOV with content of 5 µg novobiocin are designed for basic differentiation of beta-plasmacoagulase negative staphylococcus.

**PRINCIPLE**

Staphylococci are differentiated according to resistance to disc with content of novobiocin.

**DIAGNOSTIC KIT CONTAINS**

- tube with diagnostic discs - min. 500 pcs (Ref. 6025-500),  
min. 50 pcs (Ref. 6025-50)

**WARNINGS AND SECURITY PRECAUTIONS**

- **For *in vitro* diagnostics use.**
- **For professional use only**
- Follow the instructions exactly!
- All samples and inoculated products have to be considered as potentially infectious and obvious safety valid regulations must be respected while handling.
- Do not use product out of expiration.
- Before use, check if the packaging is intact. Do not use damaged kit.
- Leftovers of expired or devaluated product liquidate according to regulations of your country. Packaging materials are recyclable.

**STORAGE**

Keep antibiotic discs in dark at temperatures between +2 and +8 °C

**REQUIRED MATERIALS**

- Petri dish with Mueller-Hinton agar
- tampons, loops, burner, tubes and other basic equipment of microbiological laboratory

**RECOMMENDED PROCEDURE**

Prepare suspension of tested strain in saline (cca 0,5 McF). Spread suspension on whole area of Petri dish with sterile tampon soaked in prepared suspension. Place sterilely disc NOV on inoculated Petri dish, press it gently, so that whole surface of disc is in contact with medium. Incubate for 18-24 h at temperature 37°C.

**EVALUATION AND INTERPRETATION**

Inhibition zone < 16 mm - species of staphylococcus resistant to novobiocin.

Functionality of diagnostic discs can be verified by these control strains:

Bacterial strain	CCM / ATCC	Inhibition zone
<i>S. cohnii ssp. urealyticum</i>	4296 / 49331	< 16 mm
<i>Staphylococcus epidermidis</i>	4418 / 12228	≥ 16 mm


CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420549491430,

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