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Dátum poslednej revízie: 20. 5. 2022

**SUMMARY AND EXPLANATION**

Diagnostic discs OPT with content of 5 µg optochin are designed for preliminary identification of *Streptococcus pneumoniae* on base of its sensitivity to optochin.

PRINCIPLE

Streptococcus pneumoniae is sensitive to optochin diffused around disc, where it creates defined inhibition zone against to other α-haemolytic streptococci, which are more resistant.

DIAGNOSTIC KIT CONTAINS

- tube with diagnostic discs - min. 500 pc (Ref. 6024-500)
- min. 50 pc (Ref. 6024-50)

WARNINGS AND SECURITY PRECAUTIONS

- **For in vitro diagnostics use.**
- **For professional use only**
- Follow the instructions exactly!
- All samples and inoculated products have to be considered as potentially infectious and obvious safety valid regulations must be respected while handling.
- Do not use product out of expiration.
- Before use, check if the packaging is intact. Do not use damaged kit.
- Leftovers of expired or devaluated product liquidate according to regulations of your country. Packaging materials are recyclable.

STORAGE

Keep discs in dark at temperatures between +2 and +8 °C.

REQUIRED MATERIALS

- Petri dish with blood agar
- tampons, loops, burner, tubes and other basic equipment of microbiological laboratory

RECOMMENDED PROCEDURE

Spread suspension on whole area of Petri dish with sterile tampon soaked in suspension of tested strain in saline. Place sterilely disc OPT on inoculated Petri dish, press it gently, so that whole surface of disc is in contact with medium. Incubate for 18-24 h in CO₂ atmosphere at temperature 37°C.

EVALUATION AND INTERPRETATION

Inhibition zone ≥ 14 mm - *Streptococcus pneumoniae*


Inhibition zone < 14 mm - other α-haemolytic streptococci

Functionality of diagnostic discs can be verified by these control strains:

| Bacterial strain | CCM / ATCC | Inhibition zone |
|---------------------------------|--------------|-----------------|
| <i>Streptococcus pneumoniae</i> | 4501 / 49619 | ≥ 14 mm |
| <i>Streptococcus sp.</i> | 4685 / - | < 14 mm |

CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420549491430, e-mail: ccm@sci.muni.cz

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