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Dátum poslednej revízie: 20.05.2022

**SUMMARY AND EXPLANATION**

Diagnostic discs X-factor, V-factor a XV-factor are designed to species differentiation of *Haemophilus sp.*

**PRINCIPLE**

Paper discs are impregnated with growth factors in form of hemine (X-factor), nicotinamide adenine dinucleotide - NAD (V-factor) or with both (XV-factor). Growth factors are releasing from disc to agar. Different species of *Haemophilus* create zone around particular disc according to nutrition requirements.

**DIAGNOSTIC KIT CONTAINS**

- tube with diagnostic discs
- X-faktor (Ref. 6026-500 - min. 500ks, Ref. 6026-50 - min. 50ks)
- V-faktor (Ref. 6027-500 - min. 500ks, Ref. 6027-50 - min. 50ks)
- XV-faktor (Ref. 6028-500 - min. 500ks, Ref. 6028-50 - min. 50ks)

**WARNINGS AND SECURITY PRECAUTIONS**

- **For *in vitro* diagnostics use.**
- **For professional use only**
- Follow the instructions exactly!
- All samples and inoculated products have to be considered as potentially infectious and obvious safety valid regulations must be respected while handling.
- Do not use product out of expiration.
- Before use, check if the packaging is intact. Do not use damaged kit.
- Leftovers of expired or devaluated product liquidate according to regulations of your country. Packaging materials are recyclable.

**STORAGE**

Store diagnostic discs in dark, discs X-factor at temperature from +2 to +8 °C, discs V-factor and XV-factor at temperature from -18 to -10 °C. Used discs can be stored in refrigerator.

**REQUIRED MATERIALS**

- Petri dish with agar
- tampons, loops, burner, tubes and other basic equipment of microbiological laboratory

**RECOMENDED PROCEDURE**

Spread 24 hour culture of tested strain on agar without blood (for example Mueller Hinton agar). Place diagnostic disc with sterile loop on agar, so that whole area of disc is in contact with medium. Distance between discs must be at least 2,5 cm. Incubate 18-24 h in CO<sub>2</sub> atmosphere at temperature 35-37°C.

**EVALUATION AND INTERPRETATION**

Factor requirements. As positive result (+) is interpreted grow of colonies around diagnostic discs seen as round zone. In case of negative (-) result grow of colonies around disc is not visible.

Bacterial strain	X	V	XV
<i>Haemophilus aegyptius</i>	-	-	+
<i>Haemophilus aphrophilus</i>	-	-	+
<i>Haemophilus ducreyi</i>	+	-	+
<i>Haemophilus haemolyticus</i>	-	-	+
<i>Haemophilus influenzae</i>	-	-	+
<i>Haemophilus parahaemolyticus</i>	-	+	+
<i>Haemophilus parainfluenzae</i>	-	+	+

Functionality of diagnostic discs can be verified with this control strain:

Bacterial strain	CCM / ATCC	X + V	XV
<i>Haemophilus influenzae</i>	4457 / 49777	+	+

CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420549491430, e-mail: [ccm@sci.muni.cz](mailto:ccm@sci.muni.cz)

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