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**EN****SUMMARY AND EXPLANATION**

HIP is intended for the detection of hippurate hydrolase. The detection strip shows enzymatic activity via formation of blue color after adding HIP reagent. It is recommended for the presumptive identification of Group B streptococci, *Campylobacter jejuni* and *Gardnerella vaginalis*.

**PRINCIPLE**

Sodium hippurate as substrate of the reaction is enzymatically broken down forming a blue color, which indicates positive reaction – after HIP reagent addition.

**DIAGNOSTIC KIT CONTAINS:**

- 50 test strips for 50 determinations
- information leaflet

**WARNINGS AND PRECAUTIONS**

- **For in vitro diagnostics use.**
- **For professional use only.**
- Follow the instructions exactly!

- Used strips should be considered as potentially infectious and obvious safety valid regulations must be respected when handling.
- Observe common safety measures according to the regulations of your country.
- Before use, check if the packaging is intact. Do not use damaged kit.

**STORAGE**

At dark place at +2 to +8°C. Expiration date is placed on the label of each package.

**REQUIRED MATERIALS**

- HIP reagent (Ref. 3006)
- Tube with 0,5 ml saline

**RECOMMENDED PROCEDURE:**

1. Use pure 24 hour bacterial culture. Take a tube with 0,5 ml saline. Prepare a homogenous suspension with a density of turbidity more than 3 McF and homogenize well.
2. Insert the detection zone of strip into bacterial suspension and incubate at 35 ± 2°C for 4 - 24 hours.
3. Add 8 drops of HIP reagent into the tube.
4. Evaluate after 5-10 minutes.

**INTERPRETATION OF REACTION**

**Blue or blue-purple colour indicates positive result of the HIP test.**

**For the need of own proof function, use following recommend bacterial strains:**

Bacterial strain	CCM / ATCC	HIP	Colour
<i>Streptococcus agalactiae</i>	6187 / -	+	blue / blue-purple

CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420549491430, e-mail: [ccm@sci.muni.cz](mailto:ccm@sci.muni.cz)

TIP: use negative control – follow recommended procedure, but without bacterial culture.

Kit components doesn't contain hazardous substances.

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