## SK

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## EN

### SUMMARY AND EXPLANATION

The ONP is detection strip for detection of  $\beta$ -galactosidase. Detection strip shows colour change to yellow colour depending on metabolism of examined sample. Strip allows determination of TDA (tryptophan deaminase) and IND (indol test).

### PRINCIPLE

ONP - Breakdown of the substrate is attended by changing of colour of strip to yellow caused by releasing o-nitrophenol, this indicates positive result.

TDA - prove of tryptophan deaminase by added PHE reagent is attended by colour change - creation of red-brown colour - indicating positive result of test.

IND - Breakdown of the substrate L-tryptophan is after adding IND reagent attended by colour change to pink colour - indicating positive result.

### DIAGNOSTIC KIT CONTAINS:

- 50 test strips for 50 determinations, Information leaflet

### WARNINGS AND SAFETY PRECAUTIONS

- **For in vitro diagnostics.**
- **For professional use only**
- **Follow the instructions exactly!**
- Used strips should be considered as potentially infectious and obvious valid regulations must be respected when handling.
- Observe common safety measures according to the regulations of your country.
- Before use, check if the packaging is intact. Do not use damaged kit.

### STORAGE

At dark place at +2 to +8°C. Expiration date is placed on the label of each package.

### REQUIRED MATERIALS

Tube with 0,5 ml of saline

### RECOMMENDED PROCEDURE:

1. Use pure 24 hour bacterial culture. Prepare a homogenous suspension with density turbidity from 2 to 3 McF.
2. Insert the detection strip and dip detection zone into bacterial suspension.
3. Incubate at 35±2°C for 18-24 hours. In case of metabolically more active species, test can be evaluated after 4 hours.
4. Evaluate result of ONP test.
5. a) For proof of IND, add 2 drops of IND reagent and evaluate according to table.  
b) For proof of TDA, add 1 drops of PHE reagent and evaluate according to table.

**Do not evaluate test IND and TDA parallel in one test tube!**

### INTERPRETATION OF THE REACTION

Bacterial strain	ONP	IND	TDA
<i>Escherichia coli</i> , <i>Citrobacter koseri</i>	+ (yellow)	+ (pink)	- (yellow)
<i>Salmonella sp.</i>	- (without colour change)	- (without colour change)	- (yellow)
<i>Proteus vulgaris</i> , <i>Morgalella morgani</i> , <i>Providencia stuartii</i>	- (without colour change)	+ (pink)	+ (red-brown)
<i>Proteus penneri</i> , <i>Proteus mirabilis</i>	- (without colour change)	- (without colour change)	+ (red-brown)

### For the need of own proof function, use following recommend bacterial strains:

Control strains	CCM/ATCC	ONP	IND	TDA
<i>Escherichia coli</i>	4225 / 35218	+ (yellow)	+ (pink)	- (yellow)
<i>Proteus sp.</i>	1799 / -	- (turbidity of suspension)	+ (pink)	+ (red-brown)

CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420549491430, e-mail: [ccm@sci.muni.cz](mailto:ccm@sci.muni.cz)

TIP: use negative control – follow recommended procedure, but without bacterial culture.

Kit components doesn't contain hazardous substances.

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