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EN**SUMMARY AND EXPLANATION**

GN 24 is standardized identification system for common species of Gram negative bacteria, it is based on 24-29 miniaturized biochemical tests and internet database. List of all microorganisms, for which is kit determined, is placed at the end of leaflet.

PRINCIPLE

Kit GN 24 consists of 24 wells strip of microtitration plate in classic 96 wells format containing dehydrated substrates, where GN 24 sp is in the form of strippable microtitration plates and GN 24 fp is in the form of undivided full plate. Reconstitution of substrates is carried by inoculation of bacterial suspension. During incubation occurs colour change in wells because of metabolic activity of microorganisms. Evaluation of test results could be done by automatic reader or visually on base of colour scheme, or by colour description stated in leaflet. Results of identification can be obtained from evaluation table or by help of evaluation software, which can be found at www.diagnostics.sk/idmicro.

KIT CONTAINS - ACCORDING TO PACKAGING 40 TESTS (sp) / 100 TESTS (fp)

- 10 / 25 microtitration plates of GN 24
- 40 / 100 result forms
- 10 / 25 incubation packets
- 1 information leaflet

By interpretation of results anamnesis of a patient, source of sample, morphology of colony, microscopic morphology of strain and if necessary, results of all performed tests, mainly results of antibiogram must be considered.

REQUIRED MATERIAL**Reagents:**

- Unbuffered saline 0,85 %, 3,5-5 ml
- Paraffin oil (Ref. 3001)
- PHE reagent (Ref. 3007)
- PHS reagent (Ref. 3008)
- IND reagent (Ref. 3002)
- NIT reagent (Ref. 3005)
- VP and VP reagent (Ref. 2004 and 3004)
- Zn (Ref. 5001)
- OXI (Ref. 2001)
- PYR and PYR reagent (Ref. 2003 and 3003)
- Identification software (on website of the company)

Materials:

- Pipettes
- Tampons, loops, burner, tubes and other basic laboratory equipment.

WARNINGS AND SAFETY PRECAUTIONS

- For in vitro diagnostics use and microbial control
- For professional use only
- Follow the instructions exactly!
- Used strips should be considered as potentially infectious and this must be respected when handling.
- Observe common safety measures according to the regulations of your country.
- Before use, check if the packaging is intact. Do not use damaged kit.

STORAGE CONDITIONS

Diagnostic kits are delivered in multilayer sachets based on aluminium and organic polymers. Part of each packet is silica gel desiccant. Store kits

at temperature from +2 to +25°C. Expiration date is placed on each packaging (2 years).

Put unused microtitration strip in packed Al sachet with original silica gel desiccant, close sachet carefully and store at room temperature. Product can be stored for two weeks after opening of original Alu sachet.

SAMPLES

Isolate microorganisms, which have to be identified, from suitable unselective cultivation medium (for example blood agar, trypton - soya agar etc.) according to standard microbiological techniques.

Make Gram staining and microscopy of pure culture.

Make test of cytochromoxidase – OXI (eventually catalase test - CAT) Identify confirmed isolates on GN 24 kit. These tests are not a part of kit.

RECOMMENDED PROCEDURE

Preparation of inoculum

- Open the tube with saline or use sterile 0,85% solution of sodium chloride in deionized / distilled water.
- Take same well isolated colonies by inoculation loop from pure 18-24 hours old culture.
- Prepare well homogenized suspension of turbidity of 2 McF. Use suspension immediately after preparation.

TIP: Eventually make proof of purity by braying of same loop or tampon, by which was suspension done.

Prepared Petri dish can be used next day to make other additional tests!

Preparation of microtitration plate

- Prepare microtitration plate.
- Mark strips with numbers of examined cultures.

TIP: When using GN 24 sp, take unused strips, put it to original aluminium packet with desiccant and close carefully.

Save frame of microtitration plate after incubation for next use.

Inoculation

- Inoculate 0,1 ml of well homogenized suspension into each well of strip.
- Cover tests from URE to LYS (wells from H to C) with three drops of paraffin oil. Cover test GLU (well G) with four drops of paraffin oil.

Incubation

- Put the microtitration plate to packed PE packet, then bend end of packet under plate – this avoid drying of bacterial suspension.
- Incubate by common atmosphere at temperature 35 ± 2 °C for 18 - 24 hours.

TIP: For optimal incubation conditions, make higher humidity in incubator, for example by adding beaker of clear water or incubate in incubator with managed humidity.

In suspicion for presence of nonfermenting rods, evaluate test GLU after 4 hours. Test GLU can show positive results after 24 hours of incubation by some of metabolically more active rods. In identification program use result of GLU obtained after 4 hours of incubation.

In case of unverifiable colouring of tests in wells corresponding to carbohydrates – elongate incubation of nonfermenting bacteria (test GLU is after 4 - 24 hours negative!) to 48hours.

EVALUATION AND INTERPRETATION

Evaluate the tests after incubation with help of interpretation table, colour scheme, or by results of control strain.

Tests NAG / IND, GLR / ESL, bGL / PHE, bGA / NIT and GGT / PHS are bifunctional so after evaluation of primary reaction another result can be obtained from evaluated well of microtitration plate.

Make following bifunctional test:

Well H2: NAG / IND - add 2-3 drops of IND reagent and wait 1-2 minutes for colouring of test.

Well H3: GLR / ESL - evaluate positive test GLR for *E. coli*. Mark the negative result of the test GLR also.

Colour	Result GLR	Result ESL
No colour / creamy	-	-
Yellow	+	-
Brown	-	+

Test does not require reagent addition! Write both results GLR and ESL separately, not only result of positive test!

Well A1: bGL / PHE - add 1 drop of reagent PHE and evaluate.

Well A3: bGA / NIT - add 2 drops of NIT reagent and evaluate. In case of negative reaction of test NIT add Zinc dust to well (on tip of 1 µl loop). If it shows red colour till 10 minutes - test NIT is negative.

For nonfermenting bacteria make also following bifunctional tests:

Well A2: GGT / PHS - add 1-2 drops of PHS reagent and evaluate.

- Write down result of bifunctional tests to result form or to evaluating software.
- In case of identification of gram negative, oxidase negative, glucose nonfermenting bacteria write down to software results of test bHEM (beta hemolysis), YEP (yellow pigment production) and G42 (42°C growth), this significantly increases identification of this group of microorganisms.

IDENTIFICATION

Result of identification can be obtained by:

- identification table or identification software

Identification by identification table:

Compare results of tests and make identification by interpretation table (page 6) and by results of tests in identification table (page 7).

Identification by identification software:

Enter results of individual tests. When some of the test cannot be evaluated, it can be excluded. Software allows inputting additional tests and so it increases identification effectiveness.

Software microID is available free for our customers on the website of the company.

Identification by reader.

QUALITY CONTROL

Quality of produced diagnostic kits is systematically controlled. Chemicals are bought exclusively from certified companies and quality of these chemicals is confirmed by analytical certificate. The functionality of the kits is tested by control collection strains, controlled and tested is also present of bacterial contamination. Kits are exposed to tests of higher temperatures and samples of each batch are saved for right advisement of later reclamations.

For the need of own proof function, use recommended bacterial strains on the page 5.

CONSTRAINTS OF METHOD AND MOST OFTEN CAUSES OF WRONG IDENTIFICATION

- Diagnostic kit GN 24 is determined for identification of bacteria's named in this leaflet only.
- Only pure culture of microorganism can be used.
- Tests were not covered by paraffin oil.
- Contamination of wells by inoculum of next strip.
- Used culture is atypical strain.
- Some point of leaflet was not kept.

DETERMINATION CHARACTERISTICS

150 collection strains, strains of clinical background, and veterinary strains mentioned in database were tested:

GLU fermenting bacteria

- 92 / 88 % strains were identified correctly with additional tests
- 89 / 87 % strains were identified correctly without additional tests

GLU nonfermenting bacteria

- 96 / 87 % strains were identified correctly with additional tests
- 81 / 74 % strains were identified correctly without additional tests

Note: Data before / after slash response to genus / species identification.

WASTE LIQUIDATION

Work with material as with potentially infectious agents. Liquidate leavings according to internal procedures and legislative of your country. Kit components does not contain dangerous chemicals.

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