## SK

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**EN****SUMMARY AND EXPLANATION**

Diagnostic kit YST 8 is a standardized identification system for common identification of clinically significant species of yeast, which is based on 8 miniaturized biochemical tests. List of all microorganisms, for which is this kit determined, is placed at the end of the leaflet.

**PRINCIPLE**

The diagnostic kit YST 8 consists of 8 wells in the monostrip of the microtitration plates in classic 96 well format containing dehydrated substrates. Reconstitution of substrates runs by the inoculation of a bacterial suspension. During the incubation occurs colour change in the well because of metabolic activity of microorganisms. Evaluation of the tests results is visual on base of colour scheme or by the colour formulation described in the leaflet. Results of identification are obtained from the evaluation table with overview of profiles.

**KIT CONTAINS**

- 5 microtitration plates of YST 8
- 60 result forms
- 5 incubation bags
- 1 information leaflet

**REQUIRED MATERIAL****Reagents:**

- Paraffin oil (Ref. 3001)
- PHE reagent (Ref. 3007)

**Materials:**

- Pipettes, tampons, loops, burner, tubes and other basic laboratory equipment.

**WARNINGS AND SAFETY PRECAUTIONS**

- **For *in vitro* diagnostics use and microbial control.**
- **For professional use only.**
- Follow the instructions exactly!
- Used strips should be considered as potentially infectious and this must be respected when handling.
- Observe common safety measures according to the regulations of your country.
- Do not use after expiration.
- Before use, check if the packaging is intact. Do not use damaged kit.

By interpretation of results anamnesis of patient, source of sample, morphology of colony, microscopic morphology of batch and if necessary, results of all overrun tests, mainly results of antibiogram have to be considered.

**STORAGE CONDITIONS**

Diagnostic kits are delivered in multilayer packets made of organic polymers and aluminium. Part of each packet is silica gel desiccant. Store kits at temperatures from +2 to +25°C.

Expiration date is placed on each package.

Put unused leftover of microtitration plate in packed Al packet with original silica gel desiccant, close packet carefully and save at laboratory temperature. The product can be stored two weeks in such conditions (or to date of expiration, in case when expiration datum passes earlier).

**SAMPLES**

Isolate microorganisms, which have to be identified, from suitable unselective cultivation medium (e.g. Sabouraud agar) according to standard microbiological techniques. Consider the purity and morphology of microbial culture.

**RECCOMENDED PROCEDURE****Preparation of inoculum**

- Open the tube with whatever sterile 0,85% saline solution.
- Take same well isolated colonies by inoculation loop from 18-24 hours old culture.
- Turbidity of well homogenized suspension must be 3-5 McF.
- This suspension must be used immediately after preparation

**Inoculation**

- Mark strips with numbers of examining cultures.

- Inoculate by 0,1 ml of well homogenized suspension into each well of monostrip.
- Cover tests from URE to RAF (from H to E) with 2 – 3 drops of paraffin oil.

**Incubation**

- Put the frame with inoculated strips to packed PE bag, then bend end of the bag under the plate – this avoid dehumidifying of bacterial suspension.
- Incubate aerobically at temperature 30 ± 2 °C for 24 hours

**EVALUATION AND INTERPRETATION**

Evaluate the strip after incubation with help of evaluating form, colour scheme, or by results of control batch.

Test bGL / PHE is bifunctional so after evaluating of primary reaction another result can be obtained from same well.

**Hole A: bGL / PHE** – add 1 drop of PHE – reagent and evaluate.  
Mark results of tests to result form.

**IDENTIFICATION**

Result of identification can be obtained by:

- **Identification by identification table:**

Compare results of the tests and make identification by results of tests in identification table (page 4).

- **Identification by identification software:**

Enter results of individual tests to identification software and evaluate.

**QUALITY CONTROL**

The quality of diagnostic kits is systematically controlled. Chemicals are bought only from ISO certified companies and quality of these chemicals is confirmed by analytical certificate. The functionality of the kits is tested by collection of control strains, controlled and tested is also present of bacterial contamination. The kits are exposed to higher temperatures and samples from each batch are saved for right advisement of later reclamations. For the need of own proof function, use recommend bacterial strains (page 4).

**CONSTRAINTS OF METHOD AND MOST OFTEN CAUSES OF WRONG IDENTIFICATION**

- Some point of the leaflet was not kept
- Contamination of well by inoculum of next strip.  
Used culture is an atypical strain

**DETERMINATION CHARACTERISTICS**

Tested were 150 collection yeast strains and strains of clinical origin belonging to species included in the database

**Internal testing**

- 90 % strains were identified correctly (with or without additional tests).
- 10 % were not identified.
- 0 % was identified wrong.

**Independent testing:**

- 94 % strains were identified correctly (with or without additional tests).
- 6 % were not identified.
- 0 % was identified wrong.

**WASTE LIQUIDATION**

Work with material as with potentially infectious agents. Liquidate waste according to internal procedures and directives in accordance with the legislation of your country. Kit components does not contain dangerous chemicals.

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