## SK

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# EN

## SUMMARY AND EXPLANATION

The ENC 8 is a standardized identification system for common species of clinically significant genus *Enterococcus*. It is based on 9 miniaturized biochemical tests and internet database. List of all microorganisms, for which is this kit determined, is placed at the end of the leaflet.

## PRINCIPLE

The diagnostic kit ENC 8 consists of 8 wells in the monostrip of the microtitration plates in classic 96 well format containing dehydrated substrates. Reconstitution of substrates runs by the inoculation of a bacterial suspension. During the incubation occurs colour change in the well because of metabolic activity of microorganisms. Evaluation of the tests results is visual on base of colour scheme, or by the colour formulation described in the leaflet. Results of identification are obtained from the evaluation table with overview of profiles.

## KIT CONTAINS - 60 TESTS

- 5 microtitration plates of ENC 8
- 60 result forms
- 5 incubation bags
- 1 information leaflet

## REQUIRED MATERIAL

### Reagents:

- Paraffin oil (Ref. 3001)
- PYR (2003) and PYR reagent (3003) for pyrrolidonylarylamidase determination
- PHS reagent (Ref. 3008)

### Materials:

- pipettes, tampons, loops, burner, tubes and other basic laboratory equipment.

## WARNINGS AND SAFETY PRECAUTIONS

- **For *in vitro* diagnostics use and microbial control.**
- **For professional use only.**
- Follow the instructions exactly!
- Used strips should be considered as potentially infectious and this must be respected when handling.
- Observe common safety measures according to the regulations of your country.
- Do not use after expiration.
- Before use, check if the packaging is intact. Do not use damaged kit.

By interpretation of results anamnesis of patient, source of sample, morphology of colony, microscopic morphology of batch and if necessary, results of all overrun tests, mainly results of antibiogram have to be considered.

## STORAGE CONDITIONS

Diagnostic kits are delivered in multilayer packets made of organic polymers and aluminium. Part of each packet is silica gel desiccant. Store kits at temperatures from **+2 to +25°C**. Expiration date is placed on each package. Put unused leftover of microtitration plate in packed Al packet with original silica gel desiccant, close packet carefully and save at laboratory temperature. The product can be stored two weeks in such conditions (or to date of expiration, in case when expiration datum passes earlier).

## SAMPLES

Isolate microorganisms, which have to be identified, from suitable unselective cultivation medium (e.g. blood agar) according to standard microbiological techniques. Confirm the genus of *Enterococcus* by the agglutination, selective cultivation or by the PYR.

Confirmed isolates identify on ENC 8. Consider the purity and morphology of bacterial culture. Review haemolytic activity, motility and yellow pigment production.

## RECOMMENDED PROCEDURE

### Preparation of inoculum

- Open the tube with whatever sterile 0,85% saline solution.
- Take same well isolated colonies by inoculation loop from 18-24 hours old culture.

- Turbidity of well homogenized suspension must be **2 – 3 McF**.
- This suspension must be used immediately after preparation.

### Preparation of strip

- Prepare empty frame of microtitration plate and place there necessary number of strips.
- Mark strips with numbers of examining cultures.

### Inoculation

- Inoculate by 0,1 ml of well homogenized suspension into each well of monostrip.

- Cover test **ARG** (well H) with **2 – 3 drops of paraffin oil**.

### Incubation

- Put the frame with inoculated strips to packed PE bag, then bend end of the bag under the plate – this avoid dehumidifying of bacterial suspension.
- Incubate by common atmosphere at temperature **35 ± 2 °C** for **18 - 24 hours**.

## EVALUATION AND INTERPRETATION

Evaluate the strip after incubation with help of evaluating form, colour scheme, or by results of control batch. Test GLR / PHS is bifunctional so after evaluating of primary reaction another result can be obtained from same well.

**Well A: GLR / PHS** – add 1 drop of PHS – reagent and evaluate.

## IDENTIFICATION

Result of identification can be obtained by:

- identification table
- octale code

### Identification by identification table:

Compare results of the tests by evaluation table and make identification by results of tests in identification table (page 4).

### Identification by octale code:

Divide the tests into consecutive triples and assign to the positive test in each triple following value: first test = 1, second test = 2, third test = 4. By counting the results for each triple three-number octale code is created, which can be searched in the identification table for determination result.

## QUALITY CONTROL

The quality of diagnostic kits is systematically controlled. Chemicals are bought only from ISO certified companies and quality of these chemicals is confirmed by analytical certificate. The functionality of the kits is tested by collection of control strains, controlled and tested is also present of bacterial contamination. The kits are exposed to higher temperatures and samples from each batch are saved for right advisement of later reclamations. For the need of own proof function, use recommend bacterial strains (page 4).

## CONSTRAINTS OF METHOD AND MOST OFTEN CAUSES OF WRONG IDENTIFICATION

- The diagnostic set ENC 8 is determined for identifications of bacteria's named in this leaflet only.
- Only pure culture of microorganism can be used.
- Tests were not covered by paraffin oil.
- Contamination of wells by inoculum of next strip.
- Used culture is an atypical strain.
- Some point of the leaflet was not kept.

## DETERMINATION CHARACTERISTICS

Tested were 105 collection control strains and strains of clinical origin belonging to species included in the database:

### Internal testing

- 98 % strains were identified correctly (with or without additional tests).
- 2 % were not identified.
- 0 % was identified wrong.

### Independent testing:

- 97 % strains were identified correctly (with or without additional tests).
- 3 % were not identified.
- 0 % was identified wrong.

## WASTE LIQUIDATION

Work with material as with potentially infectious agents. Liquidate waste according to internal procedures and directives in accordance with the legislation of your country. Kit components do not contain dangerous chemicals.

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Profil after 18 - 24 hours of incubation on blood agar. Control strains serves to check functionality of individual tests only, not for proof control of identification.

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**Literatúra / Literature / Litaretura :**

1. FACKLAM, R. R., SAHM, D. F. *Enterococcus*. Manual of Clinical Microbiology, 6.vydání. American Society for Microbiology Press, Washington D.C., 1995. s. 308 – 314.
2. MANERO, A., BLANCH, A. R. Identification of *Enterococcus* spp. with a Biochemical Key. Applied and Environmental Microbiology. 1999. s. 4425 – 4430.
3. 20. Devriese L A, Collins M D, Wirth R. The genus *Enterococcus*. In: Ballows A, Trüper H G, Dworkin M, Harder W, Schleifer K H, editors; Ballows A, Trüper H G, Dworkin M, Harder W, Schleifer K H, editors. The prokaryotes. New York, N.Y.: Springer-Verlag; 1991. pp. 1465–1477.
4. Devriese L A, Pot B, Collins M D. Phenotypic identification of the genus *Enterococcus* and differentiation of phylogenetically distinct enterococcal species and species groups. J. Appl. Bacteriol. 1993;75: 399–408.
5. Facklam R R. Recognition of group D streptococcal species of human origin by biochemical and physiological tests. Appl Microbiol. 1972;23: 1131–1139.
6. Facklam R R, Collins M D. Identification of *Enterococcus* species isolated from human infections by a conventional test scheme. J Clin Microbiol. 1989;27: 731–734
7. Fertally S S, Facklam R. Comparison of physiologic tests used to identify non-beta-hemolytic aerococci, enterococci, and streptococci. J Clin Microbiol. 1987;25: 1845–1850
8. MacFaddin J F. Biochemical tests for identification of medical bacteria. Baltimore, Md: Williams and Wilkins Co.; 1980
9. Schleifer K H, Kilpper-Bälz R. Transfer of *Streptococcus faecalis* and *Streptococcus faecium* to the genus *Enterococcus* nom. rev. as *Enterococcus faecalis* comb. nov. and *Enterococcus faecium* comb. nov. Int J Syst Bacteriol. 1984;34:31–34.
10. Mabel S. Y., Facklam, R., R.,: New Test System for Identification of Aerococcus, Enterococcus, and Streptococcus Species. J. Clin. Microbiol., Oct. 1986, P. 607-611
11. Motlová, J.: Kam spěje taxonomie rodů *Streptococcus* a *Enterococcus*? I.část. Zprávy CEM (SZÚ, Praha) 2003; (7):298-305.
12. Motlová, J.: Kam spěje taxonomie rodů *Streptococcus* a *Enterococcus*? II.část. Zprávy CEM (SZÚ, Praha) 2003; 12(8):339-342.