SK

ÚČEL URČENIA A POPIS

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EN

SUMMARY AND EXPLANATION

The kit is designed to determine the minimum inhibitory concentration (MIC) of yeast. The MIC value corresponds to the lowest concentration of the antimicrobial agents, which leads to a visible reduction in the growth of the microorganism.

PRINCIPLE

The kits contain 96 well microplates containing descending gradient of dehydrated antimicrobial substrates forming a downward concentration gradient and other excipients needed to determine MIC values. Rehydration takes place in a suspension medium from a yeast suspension. The sample is visualized after incubation.

THE KIT CONTAINS

- 25 microplates and suspension media
- 25 PE bags
- working manual
- 13 result forms

REQUIRED BUT NOT SUPPLIED REAGENTS AND MATERIAL**Reagents :**

- Sterile unbuffered saline

Material:

- Pipettes, tampons, inoculation loops, burner, tubes and other basic laboratory equipment.

WARNINGS AND SAFETY PRECAUTIONS

- For *in vitro* diagnostics use only.
- For professional use only.
- Follow the instructions exactly.
- All samples and products should be considered as potentially infectious and this must be respected when handling.
- Observe common safety measures according to the regulations of your country.
- Before use, check if the packing is intact.
- Do not use damaged kit.

STORAGE CONDITIONS

Store kits at temperature from +2 to +25°C.

Expiration is placed on each packaging.

SAMPLES

Isolate microorganisms from suitable cultivation medium according to standard microbiological techniques.

RECCOMENDED PROCEDURE**Preparation of inoculum**

- Take some well isolated colonies by inoculation loop from pure 18-24 hours old culture cultivated on suitable non selective cultivation and put it into the tube with room temperature sterile saline.
- Prepare well homogenized suspension of turbidity of 0,5 McF.
- Inoculate 20 µl of freshly prepared yeast suspension into suspension media . Homogenise properly and use suspension immediately after preparation.
- If you use inoculator – please pour the yeast suspension into sterile lid of Petri dish.

PREPARATION OF MICROTITRATION PLATE

- Open aluminium sachet and mark microplate with number of examined culture
- Use microplate immediatly after opening the aluminium sachet.

INOCULATION

Inocuation by pippete – inoculate 100 µl of yeast suspension into each well of microplate.

INCUBATION

Put the microtitration plate to packed PE sachet. Bend the end of the bag under microtitration plate – to avoid drying of suspension. Incubate at temperature 35±2°C for 20±4 hours (according to EUCAST). In case of insufficient growth (well K does not show a positive reaction), prolong the incubation until the next day. Some slow-growing microorganisms (*C. glabrata*, *C. parapsilosis*, *C. guilliermondii*, *Cryptococcus sp.*, etc.) may show occasional growth after additional 12-24 hours of incubation, depending on the agar medium used, previous growth in the broth, etc. It is possible to increase the growth by the use of suspension medium enriched with bovine serum, which is not in accordance to the EUCAST methodology.

EVALUATION AND INTERPRETATION

Remove microtitration plate from PE incubation sachet and check microbial growth in control microwell (K). Repeat the test in case, if you do not see growth which indicates pink coloration. Evaluate the visual color change. The original blue color becomes pink due to the growth of yeast in each well (table below).

For the MIC value, subtract the first well, which is different in color from well K, forming a transition from pink to blue / violet!

AMB – read only **90 % inhibited microwells** as MIC value corresponding to AMB 8. Other antifungals– read min. **50% of inhibition** as MIC value in comparison with control microwell (K).

WASTE LIQUIDATION

Work with material as with potentially infectious agents. Liquidate leavings according to internal procedures and legislative of your country. The kit doesn't contain dangerous components.

QUALITY CONTROL

Control strains are on the page 3.

EVALUATION OF GROWTH

For evaluation of growth see the page 3.

Antimycotic dilution layout and MIC values (mg/l)

For antimycotic dilution layout and MIC values (mg/l) see the page 3.

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