SK

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EN

SUMMARY AND EXPLANATION

MIC stripped plates are intended for antimicrobial susceptibility testing bacteria on the basis of determination of MIC values. MIC (minimal inhibitory concentration) is the lowest concentration of antimicrobial agent, which visibly inhibits bacterial growth.

PRINCIPLE

The kit contains four stripped microtitration plates (48 tests). One test consists of one strip with 7 concentrations of antibiotics and one well for verification of growth control. Rehydration of test is realized by suspension media (4001). The results are evaluated after 16-20 hours of incubation.

THE KIT CONTAINS

- 4 microplates / 48 tests
- 4 PE bags
- Working manual
- 8 result forms
- 1 pc alu zip-log sachet

REQUIRED REAGENTS and MATERIAL

Reagents :

- Suspension media (Ref. 4001)
- Sterile saline

Material:

Pipettes, tampons, inoculation loops, burner, tubes and other basic laboratory equipment.

WARNINGS AND SAFETY PRECAUTIONS

- **For *in vitro* diagnostics use and microbial control only!**
- **For professional use only !**
- **Follow the instructions exactly!**
- Used strips should be considered as potentially infectious and this must be respected when handling.
- Observe common safety measures according to the regulations of your country.
- Before use, check if the packing is intact. Do not use damaged kit.

STORAGE CONDITIONS

Store kits at temperature from +2 to +25°C.

Expiration is placed on each packaging.

Once opened aluminium sachet can be stored in room temperature in alu zip-log sachet in dark dry place and must be used within exp. date.

SAMPLES

Isolate microorganisms, which have to be tested from suitable cultivation medium according mto standard microbiological techniques.

RECOMENDED PROCEDURE

Preparation of inoculum

- Take some well isolated colonies by inoculation loop from pure 18-24 hours old culture and put it into the tube with sterile saline.

- Prepare well homogenized suspension of turbidity of **0,5 McF**.
- Inoculate **60 µl** of freshly prepared bacterial suspension into suspension media (Ref. 4001). Use suspension immediately after preparation.

PREPARATION OF MICROTITRATION PLATE

- Open aluminium sachet and mark microplate with number of examined cultures.

INOCULATION

Inoculation by pipette – inoculate **100 µl** of bacterial suspension into each well of microplate.

INCUBATION

Put the microtitration plate to packed PE sachet. Bend the end of the bag under microtitration plate – to avoid drying of bacterial suspension. Incubate at temperature **35±1°C** for **16-20** hours (according to EUCAST).

EVALUATION AND INTERPRETATION

Remove microtitration plate from PE incubation sachet and check bacterial growth in microwell K (growth control). Repeat the test in case, that K microwell is negative (blue / purple). Evaluate change of original color and compare with control microwell (K). Concentration distribution in the strip and evaluation table are on page 3. **MIC value is the first microwell with color that differs from control microwell (K)**. Evaluate test manually or by reader.

In case, that you observe alternating of wells with positive and negative result, MIC value is where are two wells in row rated as negative (eagle phenomenon).

QUALITY CONTROL

Control strains are in table on page 4.

WASTE LIQUIDATION

Work with material as with potentially infectious agents. Liquidate leavings according to internal procedures and legislative of your country. The kit doesn't contain dangerous components.

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