SK

ÚČEL URČENIA A POPIS

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EN

SUMMARY AND EXPLANATION

The kits are intended for antimicrobial susceptibility testing of Gram positive bacteria, especially *Staphylococcus sp.* and *Streptococcus sp.* (MIC GP) and *Enterococcus sp.* (MIC ENC) on the basis of determination of MIC values. MIC (minimal inhibitory concentration) is the lowest concentration of antimicrobial agent, which visibly inhibits bacterial growth.

PRINCIPLE

The kits contain 96 well microplates containing descending gradient of dehydrated antimicrobial substrates and excipients. Rehydration of test is realized by suspension media (4001). The results are evaluated after 16-20 hours of incubation.

THE KIT CONTAINS

- 25 microplates / tests signed with red colour (MIC GP) or with black colour (MIC ENC)
- 25 PE bags
- working manual
- 13 result forms

REQUIRED REAGENTS AND MATERIAL

Reagents :

- Suspension media (Ref. 4001)
- Sterile saline

Material:

- Pipettes, tampons, inoculation loops, burner, tubes and other basic laboratory equipment.

WARNINGS AND SAFETY PRECAUTIONS

- **For *in vitro* diagnostics use and microbial control only!**
- **For professional use only !**
- **Follow the instructions exactly!**
- Used strips should be considered as potentially infectious and this must be respected when handling.
- Observe common safety measures according to the regulations of your country.
- Before use, check if the packing is intact.
- Do not use damaged kit.

STORAGE CONDITIONS

Store kits at temperature from +2 to +25°C.
Expiration is placed on each packaging.

SAMPLES

Isolate microorganisms, which have to be tested from suitable cultivation medium according to standard microbiological techniques.

RECOMMENDED PROCEDURE**Preparation of inoculum**

- Take some well isolated colonies by inoculation loop from pure 18-24 hours old culture and put it into the tube with sterile saline.
- Prepare well homogenized suspension of turbidity of 0,5 McF.
- Inoculate 60 µl of freshly prepared bacterial suspension into suspension media (Ref. 4001). Use suspension immediately after preparation.
- If you use inoculator – please pour the bacterial suspension into sterile lid.

PREPARATION OF MICROTITRATION PLATE

- Open aluminium sachet and mark microplate with number of examined cultures.

INOCULATION

Inoculation by pipette – inoculate 100 µl of bacterial suspension into each well of microplate.

Inoculation by inoculator – fill each well of microplate by 100 µl by sterile suspension media (Ref. 4001). Inoculate microplate by inoculator by following instructions: immerse sterile into bacterial suspension and apply bacterial suspension with circular movements into first half of microplate. Repeat procedure for the second half of the microplate.

INCUBATION

Put the microtitration plate to packed PE sachet. Bend the end of the bag under microtitration plate – to avoid drying of bacterial suspension. Incubate at temperature 35±1°C for 16-20 hours (according to EUCAST).

EVALUATION AND INTERPRETATION

Remove microtitration plate from PE incubation sachet and check bacterial growth in microwell K (growth control) Repeat the test in case, that K microwell is negative (blue/ purple) (table on page 3). Evaluate change of original color and compare with control microwell (K). **MIC value is the first microwell with color which differs from control microwell (K).** Evaluate test manually or by reader. ATB dilution layout and MIC values (mg/l) for *Staphylococcus sp.* / *Streptococcus sp.* is in table on page 3 (MIC GP) and on page 4 for *Enterococcus sp.* (MIC ENC).

QUALITY CONTROL

Control strains are on the page 3 (MIC GP) and page 4 (MIC ENC).

WASTE LIQUIDATION

Work with material as with potentially infectious agents. Liquidate leavings according to internal procedures and legislative of your country. The kit doesn't contain dangerous components.

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