Identifikačný systém pre Gram pozitívne koky, Identifikačný systém pre korynebaktérie, Identifikačný systém pre anaeróbne baktérie
Identification system for Gram positive bacteria, Identification system for Corynebacterium, Identification system for anaerobic bacteria
Identifikační systém pro Gram pozitivní koky, Identifikační systém pro korynebaktérie, Identifikační systém pro anaerobní bakterie
System identyfikacji ziarniaków Gram-dodatnich, system identyfikacji korynebakterii, system identyfikacji bakterii beztlenowych

SK

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EN**SUMMARY AND EXPLANATION**

Diagnostic kit GP 24 is standardized identification system for common species of Gram positive bacteria – mainly *Staphylococci*, *Streptococci* and *Enterococci*, set is based on 24-26 miniaturized biochemical tests and internet database. List of all microorganisms, for which is the kit determined, is placed at the end of leaflet.

PRINCIPLE

Kit GP 24 is made of 24 wells strips of microtitration plate in classic 96 wells format containing dehydrated substrates, where GP 24 sp is in the form of strippable microtitration plates and GP 24 fp is in the form of undivided full plate. Reconstitution of substrates runs by inoculation of bacterial suspension. During incubation occurs colour change in wells because of metabolic activity of microorganisms. Evaluation of test results could be done by automatic reader or visually on base of colour scheme, or by colour description stated in leaflet. Results of identification can be obtained from evaluation table or by help of evaluation software, which can be found at www.diagnostics.sk/idmicro.

KIT CONTAINS - ACCORDING TO PACKAGING 40 tests (sp) / 100 TESTS (fp)

10 / 25 microtitration plates of GP 24
40 / 100 result forms
10 / 25 incubation packets
1 information leaflet

REQUIRED MATERIALS**Reagents:**

Unbuffered saline 3,5-5 ml
Paraffin oil (Ref. 3001)
PHS reagent (Ref. 3008)
DMACA reagent (Ref.3009)
NIT reagent (Ref. 3005)
VP and VP reagent (Ref. 2004 and 3004)
Zn (Ref. 5001)
OXI (Ref. 2001)
PYR and PYR reagent (Ref. 2003 and 3003)
HIP reagent (Ref. 3006)
Identification software on the website of the company

Materials:

Pipettes, tampons, loops, burner, tubes and other basic laboratory equipment.

WARNINGS AND SAFETY PRECAUTIONS

- **For in vitro diagnostics use and microbial control**
- **For professional use only**
- Follow the instructions exactly!
- Used strips should be considered as potentially infectious and this must be respected when handling.
- Observe common safety measures according to the regulations of your country.
- Before use, check if the packaging is intact. Do not use damaged kit.

By interpretation of results anamnesis of patient, source of sample, morphology of colony, microscopic morphology of batch and if necessary, results of all overrun tests, mainly results of antibiogram have to be considered.

STORAGE CONDITIONS

Diagnostic kits are delivered in multilayer packets on base of aluminium and organic polymers. Part of each packet is silica gel desiccant. Store kits at temperatures from +2 to +25°C. Expiration date is placed on each packaging.

Put unused microtitration strip in packed Al sachet with original silica gel desiccant, close sachet carefully and store at room temperature. Product can be stored for two weeks after opening original Alu sachet.

SAMPLES

Isolate microorganisms, which have to be identified, from suitable unselective cultivation medium (for example blood agar, etc.) according to standard microbiological techniques. Make Gram colouring and microscope pure culture. **Make catalase test.** Confirmed isolates identify on GP 24.

RECCOMENDED PROCEDURE**Preparation of inoculum**

Open the tube with 3,5-5ml of unbuffered sterile saline. Take same well isolated colonies by inoculation loop from 18-24 hours old culture (GP 24, GP 24 COR and GP 24 BAC) or 24-48 hours old culture (GP 24 AN).

Turbidity of homogenized suspension must be 3 McF.

This suspension must be used immediately after preparation.

TIP: Eventually make proof of purity by using the same loop or tampon, by which was suspension done.

This Petri dish can be used next day to make other additive tests!

Preparation of microtitration plate

- Prepare microtitration plate
- Mark strips with numbers of examining cultures.

TIP: When using GP 24 sp first time, take other strips and put it to aluminium sachet with desiccant and close carefully. Save frame of microtitration plate for next use.

Inoculation

- Inoculate by 0,1 ml of well homogenized suspension into each well of strip.
- Cover tests URE and ARG (well H1 and H2) with 2-3 drops of paraffin oil.

GP 24 AN: we recommended determination of test IND by DMACA reagent, which is with its higher sensitivity suitable for anaerobic microorganisms. Perform the determination on the filter paper according to the instructions of DMACA reagent..

Incubation

- Put the microtitration plate to packed PE sachet, then bend end of packet under plate to avoid dehumidifying of bacterial suspension.

- **GP 24 AN, GP 24 COR and GP 24 BAC:** Incubate under ambient atmosphere at temperature 35 ± 2 °C for 18 - 24 hours.

TIP: For optimal incubation increase humidity in incubator.

- **GP 24 AN:** Incubate in anaerobic atmosphere (80% N₂, 10% H₂, 10% CO₂) at temperature 35 ± 2 °C for 18 - 24 hours.

EVALUATION AND INTERPRETATION GP 24

Make CAT test! For the CAT positive species we recommend - record the result of a NOV and COA! It will increase identification results.

For the CAT negative species we recommend - record the result of a bHEM! It will increase identification results.

Evaluate the tests after incubation with help of evaluating form, colour scheme, or by results of control batch.

Tests GLR/PHS and bGL/HIP are bifunctional.—After evaluation of primary reaction add reagent and evaluate the secondary test.

Make following bifunctional tests for G+ catalase negative bacteria:

GLR / PHS and bGL/HIP and NIT monofunctional test.

Make following bifunctional tests for G+ catalase positive bacteria:

GLR / PHS and NIT monofunctional test.

Well A2: GLR / PHS - after evaluation of primary reaction add 1-2 drops of PHS reagent and evaluate test PHS.

Well H3: NIT - add 2 drops of NIT reagent and evaluate. In case of negative reaction of test NIT after adding 1-2 drops of NIT reagent, add Zinc dust to well (on tip of 1µl loop). If it shows red colour till 10 minutes test NIT is negative. Test is monofunctional.

Well A3: bGL/HIP - after evaluation of primary reaction add 2-3 drops of HIP reagent and evaluate test HIP till 10 minutes.

- Write down result of bifunctional tests to form for evaluating results or to evaluating software.

EVALUATION AND INTERPRETATION GP 24 COR and GP24 BAC

Enter result of CAT test.

Without specifying parameter above, identification of sample will not be made.

If necessary, make following bifunctional tests: GLR/PHS, bGL/HIP and NIT monofunctional test-see above.

- Mark down result of bifunctional tests to form for evaluating results or to identification software.
- Mark result of β -hemolysis

EVALUATION AND INTERPRETATION GP 24 AN

Enter results of following tests: catalase - CAT, Gram staining - GRAM, coccus morphology - COCC and sporulation - SPOR.

Sample identification will not be made without specifying the parameters stated above.

Evaluate the tests after incubation with help of evaluating form, colour scheme, or by results of control bacterial species.

In case of decolourization of some test with content of sugar, add one drop of 0,02% unbuffered solution of bromocresol red.

Well G3: ESL – positive colour of test is most intensive after 5 minutes of exposition on air.

Make following bifunctional test

GLR / PHS and monofunctional test NIT - see above.

Write down result of bifunctional test to form for evaluating results or to evaluating software

IDENTIFICATION

Result of identification can be obtained by:

- **Identification by identification table:**
Compare results of tests and make identification by interpretation table (page 8) and by results of tests in identification table (pages 8-10).

- **Identification by identification software:**
Enter results of individual tests. When some of the test cannot be evaluated, leave the test result empty. Software allows inputting additional tests and so it increases identification effectiveness. MicroID software is for customers for free on the website of the company.

QUALITY CONTROL

Quality of diagnostic kits is systematically controlled. Chemicals are bought only from certified companies and quality of these chemicals is confirmed by analytical certificate. The functionality of the kits is tested by collection of control strains, controlled and tested is also present of bacterial contamination. Kits are exposed to tests of higher temperatures and samples of each batch are saved for right advisement of later reclamations.

Recommended bacterial strains can be used for verification of ID kit. (page 7).

CONSTRAINTS OF METHOD AND MOST OFTEN CAUSES OF WRONG IDENTIFICATION

- Diagnostic GP 24 is determined for identification of bacteria's named in this leaflet only.
- Only pure culture of microorganism can be used.
- Tests were not covered by paraffin oil.
- Contamination of wells by inoculum of next strip.
- Used culture is atypical batch.
- Some point of leaflet was not kept.

DETERMINATION CHARACTERISTICS

GP 24 : 150 typical collection bacterial strains, clinical and veterinary bacterial strains mentioned in database were tested.

97 % strains were identified correctly (with or without additional tests).

3 % were not identified or were identified wrong

GP 24 COR: 97 typical collection bacterial strains, clinical and veterinary bacterial strains mentioned in database were tested

77 % strains were identified correctly (with or without additional tests).

23 % were not identified or were identified wrong.

GP 24 AN: 50 typical collection bacterial strains, clinical and veterinary strains mentioned in database were also tested:

93 % strains were identified correctly (with or without additional tests).

7 % were not identified or were identified wrong.

WASTE LIQUIDATION

Work and dispose of the material as with potentially infectious agents in accordance within internal procedures and legislative of your country.

Kit components do not contain dangerous chemicals.

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