## SK

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## EN

### SUMMARY AND EXPLANATION

COLI is a detection strip used for screening identification of *Escherichia coli*. Detection strip shows colour change by creating indigo colour on the base of metabolism of examining samples.

### PRINCIPLE

Breakdown of substrate 4-nitrophenyl-b-D glucuronide is attended by colour change and by creating of yellow colouring of detection zone, this indicates positive result of GLR test.

Breakdown of substrate L-tryptophan is attended by colour change and by creating of pink colouring of detection zone, this indicates positive result of IND test.

### DIAGNOSTIC KIT CONTAINS:

- 50 test strips for 50 determinations  
- information leaflet

### WARNINGS AND SAFETY PRECAUTIONS

**For in vitro diagnostics use.**

- **For professional use only.**
- Follow the instructions exactly!
- Used strips should be considered as potentially infectious and obvious safety valid regulations must be respected when handling.
- Observe common safety measures according to the regulations of your country.

- Before use, check if the packaging is intact. Do not use damaged kit.

### STORAGE

At dark place at +2 to +8 °C. Expiration date is placed on the label of each package.

### REQUIRED MATERIALS

IND reagent (Ref. 3002)  
Tube with 0,5 ml of saline solution

### RECOMMENDED PROCEDURE:

1. Prepare thick bacterial suspension (≥3 McF) into the tube with saline.
2. Put in diagnostic strip COLI.
3. Incubate in thermostat by 35 ± 2°C.
4. Evaluate result of GLR test after 2 hours. Extend incubation for 4h in case of negative result.
5. Add 4 drops of IND reagent after evaluation of GLR test and evaluate result of IND test after 5 minutes.

### INTERPRETATION OF REACTION

**Positive result of tests GLR (yellow) and IND (pink) indicate *Escherichia coli*.**

**For the need of own proof function, use following recommended bacterial strains:**

Bacterial strain	CCM / ATCC	GLR	IND
<i>Escherichia coli</i>	4225 / 35218	yellow	pink

CCM: Czech Collection of Microorganisms, Masaryk University Brno, Kamenice 5, 625 00 Brno, Czech Republic, tel. +420549491430, e-mail: ccm@sci.muni.cz

TIP: Use negative control – follow recommended procedure, but without bacterial culture.

Kit components doesn't contain hazardous substances.

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