



	H	G	F	E	D	C	B	A	A'
	URE	SUC	TRE	RAF	bGA	NAG	PRO	bGL	PHE
+									
-									

Identifikačná tabuľka / Identification table / Identifikační tabulka / Azonosító táblázat

Názov / Name / Název / Név	URE	SUC	TRE	RAF	bGA	NAG	PRO	bGL	PHE	mikroskopické znaky*		
										HYP	PSH	GET
<i>Candida albicans</i>	-	+	v	-	-	+	+	-	-	-	+	+
<i>Candida glabrata</i>	-	-	+	-	-	-	-	-	-	-	-	-
<i>Candida guilliermondii</i>	-	+	(-)	+	-	-	+	+	-	-	+	-
<i>Candida kefyr</i>	-	+	-	+	+	-	-	(+)	-	-	+	-
<i>Candida krusei</i>	-	-	-	-	-	-	-	-	-	-	+	-
<i>Candida lipolytica</i>	+	-	-	-	-	-	+	v	+	-	+	-
<i>Candida lusitanae</i>	-	+	v	-	-	-	+	+	-	-	+	-
<i>Candida parapsilosis</i>	-	v	-	-	-	-	+	-	-	-	+	-
<i>Candida tropicalis</i>	-	+	+	-	-	-	-	-	-	-	+	-
<i>Cryptococcus neoformans</i>	+	-	-	-	-	v	-	+	-	-	-	-
<i>Saccharomyces cerevisiae</i>	-	+	v	+	-	-	-	-	-	-	(-)	-
<i>Trichosporon asahii</i>	+	-	-	-	-	+	-	+	(+)	+	-	-
<i>Trichosporon mucoides</i>	+	+	-	v	+	+	+	+	-	+	-	-

* mikroskopické znaky / microscopic characteristics / mikroskopické znaky / mikroszkopikus tulajdonságok

DIAGNOSTICS s. r. o., Hodská 68, 924 01 Galanta, Slovenská Republika

tel.: +421-917-742927, +421-917-742931, +421-31-7886415, fax: +421-2-69202902

IČO: 35 871 181, DIČ: 2021281394, IČ DPH: SK2021281394, www.diagnostics.sk, e-mail: info@diagnostics.sk

číslo účtu: 262 070 1536/1100 (Tatra banka), IBAN: SK84 1100 0000 0026 2070 1536, SWIFT: TATRSKBX

