

Ref. / Lot.:

Dátum: / Date: / Datum:

Kmeň č.: / Strain no.: / Kmen č.: / Törzs sz.:

Sprac.: / Processed by: / Zprac.: /
Feldolgozott:

Poznámky: / Notes: / Megjegyzések:

Rozloženie ATB / ATB dilution layout [mg/l]												
	PIP	PIT	CTX	CAZ	CPZ	CPS	CEP	MER	ERT	TGC	NET	TOB
A	128 <input type="radio"/>	128/4 <input type="radio"/>	8 <input type="radio"/>	16 <input type="radio"/>	64 <input type="radio"/>	64/32 <input type="radio"/>	16 <input type="radio"/>	16 <input type="radio"/>	2 <input type="radio"/>	8 <input type="radio"/>	16 <input type="radio"/>	8 <input type="radio"/>
B	64 <input type="radio"/>	64/4 <input type="radio"/>	4 <input type="radio"/>	8 <input type="radio"/>	32 <input type="radio"/>	32/16 <input type="radio"/>	8 <input type="radio"/>	8 <input type="radio"/>	1 <input type="radio"/>	4 <input type="radio"/>	8 <input type="radio"/>	4 <input type="radio"/>
C	32 <input type="radio"/>	32/4 <input type="radio"/>	2 <input type="radio"/>	4 <input type="radio"/>	16 <input type="radio"/>	16/8 <input type="radio"/>	4 <input type="radio"/>	4 <input type="radio"/>	0,5 <input type="radio"/>	2 <input type="radio"/>	4 <input type="radio"/>	2 <input type="radio"/>
D	16 <input type="radio"/>	16/4 <input type="radio"/>	1 <input type="radio"/>	2 <input type="radio"/>	8 <input type="radio"/>	8/4 <input type="radio"/>	2 <input type="radio"/>	2 <input type="radio"/>	0,25 <input type="radio"/>	1 <input type="radio"/>	2 <input type="radio"/>	1 <input type="radio"/>
E	8 <input type="radio"/>	8/4 <input type="radio"/>	0,5 <input type="radio"/>	1 <input type="radio"/>	4 <input type="radio"/>	4/2 <input type="radio"/>	1 <input type="radio"/>	1 <input type="radio"/>	0,125 <input type="radio"/>	0,5 <input type="radio"/>	1 <input type="radio"/>	0,5 <input type="radio"/>
F	4 <input type="radio"/>	4/4 <input type="radio"/>	0,25 <input type="radio"/>	0,5 <input type="radio"/>	2 <input type="radio"/>	2/1 <input type="radio"/>	0,5 <input type="radio"/>	0,5 <input type="radio"/>	0,064 <input type="radio"/>	0,25 <input type="radio"/>	0,5 <input type="radio"/>	0,25 <input type="radio"/>
G	2 <input type="radio"/>	2/4 <input type="radio"/>	0,125 <input type="radio"/>	0,25 <input type="radio"/>	1 <input type="radio"/>	1/0,5 <input type="radio"/>	0,25 <input type="radio"/>	0,25 <input type="radio"/>	0,032 <input type="radio"/>	0,125 <input type="radio"/>	0,25 <input type="radio"/>	0,125 <input type="radio"/>
H	1 <input type="radio"/>	1/4 <input type="radio"/>	0,064 <input type="radio"/>	0,125 <input type="radio"/>	0,5 <input type="radio"/>	0,5/0,25 <input type="radio"/>	0,125 <input type="radio"/>	0,125 <input type="radio"/>	0,016 <input type="radio"/>	0,064 <input type="radio"/>	0,125 <input type="radio"/>	K <input checked="" type="radio"/>
Hodnota MIC / MIC value:												

Ref. / Lot.:

Dátum: / Date: / Datum:

Kmeň č.: / Strain no.: / Kmen č.: / Törzs sz.:

Sprac.: / Processed by: / Zprac.: /
Feldolgozott:

Poznámky: / Notes: / Megjegyzések:

Rozloženie ATB / ATB dilution layout [mg/l]												
	PIP	PIT	CTX	CAZ	CPZ	CPS	CEP	MER	ERT	TGC	NET	TOB
A	128 <input type="radio"/>	128/4 <input type="radio"/>	8 <input type="radio"/>	16 <input type="radio"/>	64 <input type="radio"/>	64/32 <input type="radio"/>	16 <input type="radio"/>	16 <input type="radio"/>	2 <input type="radio"/>	8 <input type="radio"/>	16 <input type="radio"/>	8 <input type="radio"/>
B	64 <input type="radio"/>	64/4 <input type="radio"/>	4 <input type="radio"/>	8 <input type="radio"/>	32 <input type="radio"/>	32/16 <input type="radio"/>	8 <input type="radio"/>	8 <input type="radio"/>	1 <input type="radio"/>	4 <input type="radio"/>	8 <input type="radio"/>	4 <input type="radio"/>
C	32 <input type="radio"/>	32/4 <input type="radio"/>	2 <input type="radio"/>	4 <input type="radio"/>	16 <input type="radio"/>	16/8 <input type="radio"/>	4 <input type="radio"/>	4 <input type="radio"/>	0,5 <input type="radio"/>	2 <input type="radio"/>	4 <input type="radio"/>	2 <input type="radio"/>
D	16 <input type="radio"/>	16/4 <input type="radio"/>	1 <input type="radio"/>	2 <input type="radio"/>	8 <input type="radio"/>	8/4 <input type="radio"/>	2 <input type="radio"/>	2 <input type="radio"/>	0,25 <input type="radio"/>	1 <input type="radio"/>	2 <input type="radio"/>	1 <input type="radio"/>
E	8 <input type="radio"/>	8/4 <input type="radio"/>	0,5 <input type="radio"/>	1 <input type="radio"/>	4 <input type="radio"/>	4/2 <input type="radio"/>	1 <input type="radio"/>	1 <input type="radio"/>	0,125 <input type="radio"/>	0,5 <input type="radio"/>	1 <input type="radio"/>	0,5 <input type="radio"/>
F	4 <input type="radio"/>	4/4 <input type="radio"/>	0,25 <input type="radio"/>	0,5 <input type="radio"/>	2 <input type="radio"/>	2/1 <input type="radio"/>	0,5 <input type="radio"/>	0,5 <input type="radio"/>	0,064 <input type="radio"/>	0,25 <input type="radio"/>	0,5 <input type="radio"/>	0,25 <input type="radio"/>
G	2 <input type="radio"/>	2/4 <input type="radio"/>	0,125 <input type="radio"/>	0,25 <input type="radio"/>	1 <input type="radio"/>	1/0,5 <input type="radio"/>	0,25 <input type="radio"/>	0,25 <input type="radio"/>	0,032 <input type="radio"/>	0,125 <input type="radio"/>	0,25 <input type="radio"/>	0,125 <input type="radio"/>
H	1 <input type="radio"/>	1/4 <input type="radio"/>	0,064 <input type="radio"/>	0,125 <input type="radio"/>	0,5 <input type="radio"/>	0,5/0,25 <input type="radio"/>	0,125 <input type="radio"/>	0,125 <input type="radio"/>	0,016 <input type="radio"/>	0,064 <input type="radio"/>	0,125 <input type="radio"/>	K <input checked="" type="radio"/>
Hodnota MIC / MIC value:												