

Rozloženie ATB / ATB dilution layout [mg/l]												
	PEN	CMP	TET	AMS	ERY	CLI	CIP	GEN	T/S	VAN	TEI	NFT
A	4 <input type="radio"/>	32 <input type="radio"/>	8 <input type="radio"/>	16 <input type="radio"/>	8 <input type="radio"/>	16 <input type="radio"/>	8 <input type="radio"/>	4 <input type="radio"/>	4/76 <input type="radio"/>	16 <input type="radio"/>	16 <input type="radio"/>	128 <input type="radio"/>
B	2 <input type="radio"/>	16 <input type="radio"/>	4 <input type="radio"/>	8 <input type="radio"/>	4 <input type="radio"/>	8 <input type="radio"/>	4 <input type="radio"/>	8 <input type="radio"/>	2/38 <input type="radio"/>	8 <input type="radio"/>	8 <input type="radio"/>	64 <input type="radio"/>
C	1 <input type="radio"/>	8 <input type="radio"/>	2 <input type="radio"/>	4 <input type="radio"/>	2 <input type="radio"/>	4 <input type="radio"/>	2 <input type="radio"/>	4 <input type="radio"/>	1/19 <input type="radio"/>	4 <input type="radio"/>	4 <input type="radio"/>	32 <input type="radio"/>
D	0,5 <input type="radio"/>	4 <input type="radio"/>	1 <input type="radio"/>	2 <input type="radio"/>	1 <input type="radio"/>	2 <input type="radio"/>	1 <input type="radio"/>	2 <input type="radio"/>	0,5/9,5 <input type="radio"/>	2 <input type="radio"/>	2 <input type="radio"/>	16 <input type="radio"/>
E	0,25 <input type="radio"/>	2 <input type="radio"/>	0,5 <input type="radio"/>	1 <input type="radio"/>	0,5 <input type="radio"/>	1 <input type="radio"/>	0,5 <input type="radio"/>	1 <input type="radio"/>	0,25/4,75 <input type="radio"/>	1 <input type="radio"/>	1 <input type="radio"/>	8 <input type="radio"/>
F	0,125 <input type="radio"/>	1 <input type="radio"/>	0,25 <input type="radio"/>	0,5 <input type="radio"/>	0,25 <input type="radio"/>	0,5 <input type="radio"/>	0,25 <input type="radio"/>	0,5 <input type="radio"/>	0,125/2,375 <input type="radio"/>	0,5 <input type="radio"/>	0,5 <input type="radio"/>	4 <input type="radio"/>
G	0,064 <input type="radio"/>	0,5 <input type="radio"/>	0,125 <input type="radio"/>	0,25 <input type="radio"/>	0,125 <input type="radio"/>	0,25 <input type="radio"/>	0,125 <input type="radio"/>	0,25 <input type="radio"/>	0,064/1,187 <input type="radio"/>	0,25 <input type="radio"/>	0,25 <input type="radio"/>	2 <input type="radio"/>
H	0,032 <input type="radio"/>	0,25 <input type="radio"/>	0,064 <input type="radio"/>	0,125 <input type="radio"/>	0,064 <input type="radio"/>	0,125 <input type="radio"/>	0,064 <input type="radio"/>	0,125 <input type="radio"/>	0,032/0,594 <input type="radio"/>	0,125 <input type="radio"/>	0,125 <input type="radio"/>	K <input checked="" type="radio"/>
Hodnota MIC / MIC value:												

Rozloženie ATB / ATB dilution layout [mg/l]												
	PEN	CMP	TET	AMS	ERY	CLI	CIP	GEN	T/S	VAN	TEI	NFT
A	4 <input type="radio"/>	32 <input type="radio"/>	8 <input type="radio"/>	16 <input type="radio"/>	8 <input type="radio"/>	16 <input type="radio"/>	8 <input type="radio"/>	4 <input type="radio"/>	4/76 <input type="radio"/>	16 <input type="radio"/>	16 <input type="radio"/>	128 <input type="radio"/>
B	2 <input type="radio"/>	16 <input type="radio"/>	4 <input type="radio"/>	8 <input type="radio"/>	4 <input type="radio"/>	8 <input type="radio"/>	4 <input type="radio"/>	8 <input type="radio"/>	2/38 <input type="radio"/>	8 <input type="radio"/>	8 <input type="radio"/>	64 <input type="radio"/>
C	1 <input type="radio"/>	8 <input type="radio"/>	2 <input type="radio"/>	4 <input type="radio"/>	2 <input type="radio"/>	4 <input type="radio"/>	2 <input type="radio"/>	4 <input type="radio"/>	1/19 <input type="radio"/>	4 <input type="radio"/>	4 <input type="radio"/>	32 <input type="radio"/>
D	0,5 <input type="radio"/>	4 <input type="radio"/>	1 <input type="radio"/>	2 <input type="radio"/>	1 <input type="radio"/>	2 <input type="radio"/>	1 <input type="radio"/>	2 <input type="radio"/>	0,5/9,5 <input type="radio"/>	2 <input type="radio"/>	2 <input type="radio"/>	16 <input type="radio"/>
E	0,25 <input type="radio"/>	2 <input type="radio"/>	0,5 <input type="radio"/>	1 <input type="radio"/>	0,5 <input type="radio"/>	1 <input type="radio"/>	0,5 <input type="radio"/>	1 <input type="radio"/>	0,25/4,75 <input type="radio"/>	1 <input type="radio"/>	1 <input type="radio"/>	8 <input type="radio"/>
F	0,125 <input type="radio"/>	1 <input type="radio"/>	0,25 <input type="radio"/>	0,5 <input type="radio"/>	0,25 <input type="radio"/>	0,5 <input type="radio"/>	0,25 <input type="radio"/>	0,5 <input type="radio"/>	0,125/2,375 <input type="radio"/>	0,5 <input type="radio"/>	0,5 <input type="radio"/>	4 <input type="radio"/>
G	0,064 <input type="radio"/>	0,5 <input type="radio"/>	0,125 <input type="radio"/>	0,25 <input type="radio"/>	0,125 <input type="radio"/>	0,25 <input type="radio"/>	0,125 <input type="radio"/>	0,25 <input type="radio"/>	0,064/1,187 <input type="radio"/>	0,25 <input type="radio"/>	0,25 <input type="radio"/>	2 <input type="radio"/>
H	0,032 <input type="radio"/>	0,25 <input type="radio"/>	0,064 <input type="radio"/>	0,125 <input type="radio"/>	0,064 <input type="radio"/>	0,125 <input type="radio"/>	0,064 <input type="radio"/>	0,125 <input type="radio"/>	0,032/0,594 <input type="radio"/>	0,125 <input type="radio"/>	0,125 <input type="radio"/>	K <input checked="" type="radio"/>
Hodnota MIC / MIC value:												